

Virus Transportation Media Solution

PT Elo Karsa Utama menyediakan solusi untuk pengambilan dan pengiriman spesimen virus.

VTM dapat menjaga keseimbangan osmotik & pH hasil swab sehingga spesimen tidak rusak saat diuji PCR. Untuk list media VTM terlampir

NUNC Cryotube sebagai penyimpanan spesimen dan media harus bebas dari DNase & RNase dengan tingkat sterilitas yang tinggi sehingga meningkatkan keakuratan hasil uji PCR.

Produk BioBottle dan NUNC Cryotube dirancang sesuai dengan regulasi IATA agar dapat menjaga integritas spesimen dalam suhu -20°C selama 2 hari (48 jam) dan tekanan rendah sehingga dapat memberikan keamanan saat dipindahkan melalui udara.


FEATURE

NUNC Cryotube 1.8 mL

- Bebas DNA-ase & RNA-ase, non sitotoksik dan non mutagenik material
- Sterile (SAL 10⁶) dan DNA/RNA, protein low binding Surface
- Sesuai regulasi IATA untuk pengiriman udara

BioBottle

- Durable dan memiliki sertifikasi IATA



Source: Pedoman Pencegahan dan pengendalian Coronavirus disease (Covid-19), 20 Maret 2020

Catalogue	Description
375418	Nunc Cryotube 1.8 mL, External thread, Self standing
GC2	Biobottle - BioFreeze complete
SLGP033RS	Millex Syringe Filter Unit, 0.22 µm, polyethersulfone, 33 mm
SLGV033RS	Millex Syringe Filter Unit, 0.22 µm, PVDF, 33 mm
566-0020	Nalgene Rapid Flow 0.2 µm, PES, 75 mm diam, 500 mL
169-0045	Nalgene Rapid Flow 0.45 µm, PES, 90 mm diam, 500 mL
H9269-500ML	Hanks' Balanced Salt solution (Magnesium, Calcium and Phenol red)
S8761-100ML	Sodium bicarbonate solution (7.5%)
A8412-100ML	Bovine Serum Albumin solution 7.5% in DPBS
A2942-20ML	Amphotericin B solution 250 µg/mL in deionized water
12657-5GMCN	Albumin, Bovine Serum, Fraction V, Crystalline
G1397-10ML	Gentamicin solution 50 mg/mL in deionized water
M5650-500ML	Eagle Minimum Essential Medium (EMEM) with amino acid without L-Glutamine
M0325-500ML	Eagle Minimum Essential Medium (EMEM) with amino acid and L-Glutamine
H0887-100ML	HEPES Solution
P4333-100ML	Penicillin-Streptomycin with 10,000 units penicillin and 10 mg streptomycin/mL
1.09141.1000	NaOH 1 N
M7145-100ML	Non Essential Amino Acid (100X)
F2442-500ML	Foetal Bovine Serum (US Origin)
F7524-50ML	Foetal Bovine Serum (non US Origin)



Head Office:
 Jl. Raya Kebayoran Lama 34E
 Jakarta Selatan – 12220
 Phone: +62 – 21 – 739 2856, 720 1893
 Fax: +62 – 21 – 726 0177
 Email: info@elokarsa.com

Bandung
 Istana Pasteur Regency Kav CRB-87
 Jl. Sukaraja Pasteur, Bandung
 Phone: +62 – 22 – 860 65333
 Fax: +62 – 22 – 860 65393
 E-mail: salesbdg@elokarsa.com

Yogyakarta
 Jl. Monumen Yogya Kembali No.116
 Mlati, Sleman, Yogyakarta
 Phone: +62 – 274 – 502 2016, 502 2013
 Fax: +62 – 274 – 502 2013
 E-mail: salesygy@elokarsa.com

Bogor
 Jl. Ring Road Bogor Utara No.178
 Taman Yasmin Sektor VI, Bogor
 Phone: +62 – 251 – 754 3140
 Fax: +62 – 251 – 754 3108
 E-mail: salesbgr@elokarsa.com

Surabaya
 Margorejo Indah VIII, Blok B-721,
 Surabaya - 60238
 Phone: +62 – 31 – 848 0326
 Fax: +62 – 31 – 848 0327
 E-mail: salessby@elokarsa.com