

# Extraction & Purification Kit



## Absolutely RNA Kit

- Generate highly pure RNA that is DNA-free
- Streamlined 30-minute methods save time
- Application-based testing ensures highest level of quality control
- Filter column based

Item	Amount	Catalog no.	Note
Absolutely RNA Miniprep Kit	50 preps	400800	Tissue, culture cell, nasopharyngeal samples with viral RNA, 100-600 $\mu$ L
Absolutely RNA Microprep Kit	50 preps	400805	Tissue, culture cell, nasopharyngeal samples with viral RNA, 150 $\mu$ L
Absolutely RNA 96 Microprep Kit	2 plates	400793	Tissue, culture cell, nasopharyngeal samples with viral RNA, 150 $\mu$ L
Absolutely RNA Nanoprep Kit	50 preps	400753	Tissue, culture cell, nasopharyngeal samples with viral RNA, 100 $\mu$ L
Absolutely mRNA Purification Kit	10 preps	400806	Magnetic Bead-Based
Absolutely RNA FFPE Kit	50 preps	400809	FFPE sample
Absolutely RNA miRNA Kit	50 preps	400814	Tissue and culture cell

## StrataPrep Plasmid Miniprep Kit

- Accommodates DNA from 2.9 Kb plasmids to 45 Kb cosmids
- 50 preps (Cat No. 400761) & 250 preps (Cat No. 400763)

## StrataPrep DNA Gel Extraction

- Allow rapid, high-yield isolation of pure DNA from agarose gels
- 50 preps (Cat No. 400766) & 250 preps (Cat No. 400768)

## StrataPrep PCR Purification Kit

- Accommodates purification of PCR products, from 100 bp to long PCR amplicons
- Ideal purification for high throughput sequencing and microarrays
- 50 preps (Cat No. 400771) & 250 preps (Cat No. 400773)